

# Product Data Sheet

Catalogue Number: SKU0012  
Revision: 1

## LunaGel™ Photocrosslinkable Extracellular Matrix Porcine Skin Gelatin, High Stiffness (SKU0012)

For Research Use Only. Not for use in human or animal therapeutic or diagnostic procedures.

### Product Description

LunaGel™ is a ready-to-use photocrosslinkable extracellular matrix (ECM) based on chemically modified pharmaceutical grade gelatin. The major components of LunaGel™ include ECM proteins collagen type I, II, IV, and V, as well as musculoskeletal tissue glycoproteins and proteoglycans. LunaGel™ retains the intrinsic cell-instructive bioactivity of natural ECMs, facilitating cell attachment, proliferation, differentiation, migration, and proteolytic degradation.



LunaGel™'s unique photocrosslinking technology allows unprecedented control over matrix porosity and stiffness, allowing researchers to replicate the physicochemical properties of a variety of healthy and diseased tissues in a simple 3D cell culture format. Porcine Skin Gelatin creates optically transparent hydrogels which are stable at room temperature and compatible with standard imaging systems and bioassays. Viable cells, organoids, and spheroids cultured in gelatin based LunaGel™ ECMs can be easily harvested using a LunaGel™ Cell Recovery Kit.

### Product Specification

Option	High stiffness (0 - 25 kPa), 7.5ml
Kit Contents	5 ml LunaGel™ ECM solution supplied as a sterile 1.5X stock solution in PBS 5 vials of lyophilized photoinitiator
Use	3D cell culture or musculoskeletal cell types; tissue engineering; bioprinting
Formulation	Contains ECM proteins collagen type I, III, IV, and V, as well as connective tissue glycoproteins and proteoglycans. No active growth factors present.
Physical State	Supplied as solution
pH	6.5 - 7.5
Cell Recovery	Use LunaGel™ Cell Recovery Kit
Storage	Stored at 4 - 8 °C, protected from light. Ships at ambient temperature.
Expiry	12 months from the manufacture date Following reconstitution in buffer, store the photoinitiator solution at 4 - 8 °C protected from light, and use within 7 days.

### 3D Cell Culture Workflow

LunaGel™ uses visible light polymerization to create 3D cell culture models with physiological stiffness within just minutes. Cells are resuspended in the LunaGel™ precursor solutions, dispensed into standard well plates, and cured to form stable hydrogels by exposure to cell-friendly visible light in the LunaCrosslinker™.



### Related Products

LunaCrosslinker™ (SKU 0004) – Visible Light Photocrosslinking Device  
LunaGel™ Porcine Skin Gelatin Low Stiffness Kit (SKU 0002)  
LunaGel™ Cell Recovery Kit (SKU 0015)

### Links

[LunaGel™ Safety Data Sheet \(SDS\)](#)  
[Full protocol and video demonstrations](#)

## 3D Cell Culture Protocol

### General Notes

The LunaGel™ ECM based on mammalian gelatin (bovine or porcine) may undergo reversible thermal gelation at temperatures below ~30°C. Heating the ECM solution to 37°C will liquify the solution in preparation for cell encapsulation.

The mechanical properties of LunaGel™ ECM vary depending on the type of tissue culture plastic and the volume used. We recommend optimizing the mechanical properties to suit your cell type. Try 2-, 4-,

and 8-min exposure, respectively, and observe cell growth under a microscope after 5 - 7 days.

The optimal cell concentration is cell-type dependent. As a starting point, we recommend 100,000 - 500,000 cells per ml of LunaGel™ ECM for spheroid forming cell types, 1 - 5 million cells per ml for stromal and mesenchymal cell types, and 6-8 million cells per ml for endothelial cell tube formation assays.

# Product Data Sheet

Catalogue Number: SKU0012  
Revision: 1

## Required materials and devices

- LunaCrosslinker™ - Visible Light Photocrosslinking Device
- Phosphate-buffered saline (PBS), pH 7.4
- Non-tissue culture treated polystyrene plates

## Experimental procedure for a final volume of 750 µl of cell-laden ECM\*

\*Adjust volumes as required. Preparing lower volume < 500 µl makes handling more challenging.

1. Place the 1.5X LunaGel™ ECM solution into a water bath at 37 °C for approximately 15 min or until liquid.
2. Reconstitute one vial of LunaGel™ Photoinitiator with 500 µl PBS and store protected from light.
3. Lift and count cells according to your standard protocol. Inhibit and remove trypsin/protease solution from the cells.
4. Transfer the required number of cells into a fresh reaction tube and pellet by centrifugation. Remove the entire supernatant, taking care to minimize liquid residues that may dilute the ECM solution in later steps.
5. Remove the LunaGel™ ECM solution from the water bath and mix thoroughly by pipetting up and down. Take care not to introduce air bubbles, LunaGel™ High Stiffness may be somewhat viscous – consider using a positive displacement pipette for most accurate results.
6. Add 500 µl of the LunaGel™ ECM solution to the cell pellet and gently pipette up and down to resuspend cells.
7. Add 250 µl of the photoinitiator solution to the cell suspension. Store the remaining photoinitiator solution at 4 – 8 °C, protected from light, for future use.
8. Mix thoroughly by pipetting up and down to ensure a homogenous cell suspension. Take care to avoid the introduction of air bubbles.
9. Plate the mixture in a culture dish of your choice. You can use multiwell plates or glass-bottom cell culture chambers (ideal for imaging). We recommend using non-treated culture plasticware to minimize cell adherence to tissue culture plastic. The recommended hydrogel volume for well plates is listed below.

**Pro tip:** Use reverse pipetting to avoid introducing air bubbles when dispensing the mixture into well plates.

Multiwell type	6 well plate	12 well plate	24 well plate	48 well plate	96 well plate	384 well plate
Volume per well	1,200 µl	600 µl	300 µl	150 µl	50 µl	10 µl

10. Crosslink the cell-laden LunaGel™ ECM by light exposure using the Luna Crosslinker™.
11. Add sufficient cell culture medium to cover the gel and incubate in a tissue culture incubator. The recommended cell culture media volumes are listed below.

Multiwell type	6 well plate	12 well plate	24 well plate	48 well plate	96 well plate	384 well plate
Volume per well	1,200 µl	600 µl	300 µl	150 µl	50 µl	10 µl

12. Change the cell culture medium as required, taking care not to damage the gel samples.

## Troubleshooting Guide

Problem	Solution
LunaGel™ ECM solution solidifies during protocol	LunaGel™ ECM solution solidifies during protocol
Air bubbles in LunaGel™ ECM solution.	Centrifuge solution (with or without cells) at 300g for 1 min and mix by pipetting up and down.
The LunaGel™ ECM solution does not crosslink when exposed to light.	Ensure photoinitiator solution is prepared fresh and/or extend crosslinking time.
Cells are not viable.	Reduce crosslinking time.
LunaGel™ ECM samples dissolve following cell encapsulation.	Traces of trypsin or other proteases may degrade the hydrogel samples. Ensure complete removal of trypsin before cell encapsulation by washing the cells pellet with medium of buffer.

Please contact us at [info@gelomics.com](mailto:info@gelomics.com) for questions or more information.